

Exceptional service in the national interest



Applications of high-throughput multibeam SEM in materials science

J. R. Michael, Craig Nakakura

Sandia National Laboratories, Albuquerque, NM

Tomasz Garbowski, Anna Lena Eberle, Thomas Kemen and Dirk Zeidler
Carl Zeiss Microscopy GmbH, Oberkochen, Germany




Sandia National Laboratories is a multi-program laboratory managed and operated by Sandia Corporation, a wholly owned subsidiary of Lockheed Martin Corporation, for the U.S. Department of Energy's National Nuclear Security Administration under contract DE-AC04-94AL85000. SAND NO. 2011-XXXXP

A new way to look at images! (HD View from Microsoft as an example)



Buffering...




<http://research.microsoft.com/en-us/um/redmond/groups/ivm/HDView/HDabout.htm>

How can large area high-resolution imaging be accomplished in the SEM?

Need huge numbers of pixels with small spacing between them for high resolution to cover the large area image.

This is slow!




$$T_f > \frac{(4 \times 10^{-18}) n_{pixels}}{\delta(DQE) C^2 I_B}$$

Time required to image 1 mm² with 0.5 nA beam current with 4 nm pixels.

Not including stage moves or other overhead.



Properly engineered multibeam SEM can improve these times.

Multi-beam SEM electron optics



Electron beam energy is at 30 kV throughout the probe producing column, sample is biased to result in a landing energy between 100 eV and 6 kV

Proc. SPIE, doi:
[10.1117/12.2069119](https://doi.org/10.1117/12.2069119)



High throughput data acquisition with a multi-beam SEM

Anna Lena Keller^a, Dirk Zeidler^a, Thomas Kemen^a
^aCarl Zeiss Microscopy GmbH, Carl-Zeiss-Strasse 22, 73447 Oberkochen, Germany

ABSTRACT


Conventional scanning electron microscopes are limited in their ultimate data acquisition rate at a given resolution by statistical electron-electron interaction (so-called Coulomb interaction) as well as band width of detectors and deflection systems. We increased imaging speed dramatically by using multiple electron beams in a single column and parallel detection of the secondary electrons. The multi-beam SEM generates multiple overlapping images during a single scan pass, thereby covering a larger area in shorter time as compared to a single-beam SEM at the same pixel size. This addresses the upcoming need for high speed imaging at electron microscopic resolution to investigate larger and larger areas and volumes.

Keywords: Multi-beam, SEM, high speed imaging, beam splitter

1. INTRODUCTION

Electron Source and projective column

Projective
column collects
61 separate
secondary
electron signals



Schottky electron source

Multi-optic forms 61
beams



Magnetic prism

Imaging lens

Sample


Diagram courtesy Carl Zeiss Microscopy

Zeiss MultiSEM



Zeiss MultiSEM installation completed in June, 2015


Relationship between scanned beams and secondary electron signals




Each of the 61 primary beams of electrons creates a beam of secondary electrons that are collected in parallel. Maximum throughput is 1.2 billion pixels/sec.

For large area mapping the stage is moved the process repeated to cover the desired area


MultiSEM Performance



Spot size




Beam current




Calibrations and acceptance criteria are done for each beam in the 61 beam array!

These are just a few of the many performance metrics developed to deal with multiple beams in one SEM.


Crosstalk – Occurs when secondary electron signal from one beam extends into an adjacent beam



Proper alignment of the secondary electron beams with the optical fibers reduces crosstalk to acceptably low levels




Metroboost calibration chip at 6 kV



Ghost image

Large area imaging in the multiSEM



Sample requirements:


Sample must be relatively flat – Rough or samples that are not flat can cause difficulties with the collection of the secondary electrons. Working distance of 1.4 mm.

Auto focus at each stage position or interpolation improves imaging of rough samples

Samples must be conductive – Required to get the expected landing energy, conductive coatings are helpful

Some samples that have areas that are conductive surrounded by regions of non-conductive material can be imaged without coating.


Metroboost image test chip (1.5 keV landing energy)



2.3 x 2.1 mm area acquired
at 3.8 nm/pixel.


480 mFOVs collected in 25
minutes.

Calibration sample (1.5 keV landing energy)




Metroboost calibration sample

Calibration sample (1.5 keV landing energy)



Metroboost calibration sample

Calibration sample (1.5 keV landing energy)






Single image tile

3340x2868 pixels


3.8 nm/ pixwl

SRAM test Structure (3.0 kV landing energy)



2 X 1 mm area (210 mFOVs) scanned at 3.8 nm/pixel in 32 minutes






Polished sandstone sample - coated with AuPd (1.5keV)




3 mm X 3 mm area at
6 nm/pixel

1569 mFOVs scanned
in 19 minutes

Polished sandstone sample - coated with AuPd (1.5keV)




Tin-plated copper – tin whisker mitigation study (1.5 kV Landing energy)




3 X 2.5 mm area at
3.8 nm/pixel




815 mFOVs acquired
in 51 minutes



Tin-plated copper – tin whisker mitigation study









1207 mFOVs collected in 84 minutes




1.5 keV
landing energy

Tin-plated copper – tin whisker mitigation study

Various whisker images extracted from large area




Low landing energy capabilities




6 eV landing energy

Low landing energy capabilities produce interesting imaging modes and possibilities


SRAM test structure



7 eV landing energy




100 eV landing energy




500 eV landing energy

Low voltage imaging with MultiSEM

Tin whisker



1.5 kV landing energy



2 eV landing energy

Conclusions

The need for speed:

- ✓ 61 beam multSEM – 1.22 Gpixels/sec maximum throughput
- ✓ 4 nm pixel size
- ✓ Large areas can be imaged quickly with good contrast
- ✓ Secondary electron imaging only – no BSE, x-ray etc.
- ✓ Landing energy can be varied over a useful range (few eV to 6 kV)
- ✓ Samples must be conductive and reasonably flat
- ✓ Operation similar to SEM with additional complications due to multibeams
- ✓ Interesting applications in microelectronics, geology and materials science